

Rowan University

Rowan Digital Works

Stratford Campus Research Day

25th Annual Research Day

May 6th, 12:00 AM

Safety of Silver Oxide Coated Biomaterials in Mice

Michael Klug

Rowan University

Darien L. Seidman

Rowan University

Rahyan Mahmoud

Rowan University

Destiny Morot

Rowan University

Lei Yu

Rowan University

See next page for additional authors

Follow this and additional works at: https://rdw.rowan.edu/stratford_research_day



Part of the [Animal Experimentation and Research Commons](#), [Biochemistry Commons](#), [Biomedical and Dental Materials Commons](#), [Biotechnology Commons](#), [Laboratory and Basic Science Research Commons](#), and the [Molecular Biology Commons](#)

Let us know how access to this document benefits you - share your thoughts on our [feedback form](#).

Klug, Michael; Seidman, Darien L.; Mahmoud, Rahyan; Morot, Destiny; Yu, Lei; Hettinger, Jeffrey D; and Demarest, Renee M, "Safety of Silver Oxide Coated Biomaterials in Mice" (2021). *Stratford Campus Research Day*. 68.

https://rdw.rowan.edu/stratford_research_day/2021/may6/68

This Poster is brought to you for free and open access by the Conferences, Events, and Symposia at Rowan Digital Works. It has been accepted for inclusion in Stratford Campus Research Day by an authorized administrator of Rowan Digital Works.

Author(s)

Michael Klug, Darien L. Seidman, Rahyan Mahmoud, Destiny Morot, Lei Yu, Jeffrey D Hettinger, and Renee M Demarest

Safety of Silver-Oxide Coated Biomaterials in Mice

Michael Klug⁴, Darien L. Seidman⁴, Rahyan Mahmoud³, Destiny Morot¹, Lei Yu⁵, Jeffrey Hettinger⁶
and Renee M. Demarest, Ph.D.^{1, 2, 3, 4}

Department of Molecular Biology¹, Department of Cell Biology & Neuroscience²

Rowan University Graduate School of Biomedical Sciences³, Rowan University School of Osteopathic Medicine⁴,

Department of Chemistry and Biochemistry⁵, Rowan University, Department of Physics and Astronomy, Rowan University⁶

Abstract

It has been demonstrated that silver-oxide coatings designed by our collaborators are able to prevent *E. Coli* and *P. Aeruginosa* attachment to biomaterials *in vivo*. These findings demonstrate that such coatings show promise in preventing the development of biofilm on biodevices. However, it is unknown if the use of silver-oxide in this fashion is toxic *in vivo*. The goal of this project was to determine whether our silver-oxide coatings are safe to use *in vivo*. To assess the toxicity of our silver-oxide formula, mice were implanted with either silver-oxide coated titanium discs or uncoated titanium discs. Blood samples were drawn at pre-determined time points in order to determine AST and ALT levels via ELISA assay. Preliminary results demonstrate no acute liver injury after 3 months with the discs. However, it appears as if silver is accumulating in tissues over time. Histological analysis at one year shows evidence of simple steatosis in livers. While our group is continuing to investigate the safety and efficacy of these silver-oxides coatings *in vivo*, preliminary data shows that they may have some toxicity and the silver-oxide formula may need to be altered and retested.

Background

Silver is a material that has been used in medicine for hundreds of years to promote wound healing. Silver oxidizes when exposed to fluid and its ions act as antimicrobial agents. Because of the recent increase in healthcare acquired antibiotic resistant bacterial strains^{1,2}, the use of silver coatings for medical devices is being explored extensively. Because of its low solubility, certain components in the body can act as “sponges” for silver ions, with ionization in the local microenvironment being ideal³. Biofilm formation on implanted materials or biodevices is a common issue causing severe health consequences in patients including revision surgery¹. This biofilm accumulation may be able to be combated through silver antimicrobial capabilities.

Dr. Demarest and her colleagues have developed a sputtering technique that can coat most solid materials with silver-oxide and control the elution rate by altering the sputtering formula⁴. These synthesized coatings also have three-orders-of-magnitude greater silver ion elution rates when compared to high-surface-area silver nanoparticles. Our collaborators have developed silver-oxide coatings that can prevent *E. coli*, *S. Aureus*, and *P. Aeruginosa* growth *in vitro*⁴. This silver oxide-coating is also effective against *E. Coli* and *P. Aeruginosa* attachment *in vivo*⁵. These findings demonstrate that silver-oxide coatings show initial promise in preventing the development of biofilm on biodevices. However, it is unknown if use of silver oxide in this fashion is toxic *in vivo*. Thus, the focus of this project was to determine if the silver oxide coatings are safe to use *in vivo*.

Results

In order to determine whether liver damage is occurring when mice are implanted with silver-oxide coated discs, AST and ALT Elisa assays were performed on serum. As a positive control for the assay, acetaminophen was administered to control mice in order to increase AST and ALT levels (Fig. 1). After 1 wk and 3 months post-insertion of discs the AST and ALT levels were not significantly increased in mice with silver-oxide coated discs compared to mice with uncoated discs (Fig 2). ICP-MS analysis of tissue at 1 wk post-insertion shows significant silver accumulation in some tissues and an increasing trend in many other tissues (Fig 3). Histological analysis demonstrates simple steatosis in the livers of mice with silver-oxide coated discs inserted compared to control mice with uncoated discs inserted (Fig 4).




Fig. 1. Trial ALT assay. A trial ALT was performed at 1 wk using an acetaminophen treated mice as a control in triplicate. Average values of the triplicate samples were graphed. ANOVA analysis was performed in order to determine significance. *p-value < 0.05

Fig. 2. ALT levels do not increase in mice 1-wk after disc insertion. AST and ALT were performed on serum for mice inserted with an uncoated or silver-oxide coated disc 1-wk post-insertion. ANOVA analysis was performed in order to determine significance. *p-value < 0.05; N.S. = no significant change




Fig. 3. ICP-MS analysis to determine silver accumulation in tissue. Mice were euthanized the indicated tissue was harvested 1-wk after insertion of coated or uncoated discs. Tissue was boiled in nitric acid and analyzed by ICP-MS to determine silver ppm/mg.




Fig. 4. Histological analysis to determine evidence of tissue injury. Mice were euthanized the indicated tissue was harvested 1-yr after insertion of coated or uncoated discs. Tissue was fixed in 10% neutral buffered formalin, processed, embedded in paraffin and sectioned at 5 microns. Sections were mounted on slides and H&E staining performed. Images were acquired using 10X and 40X objectives.

Materials & Methods

Disc Implantation

All protocols were approved by RowanSOM IACUC. The mice were placed in the induction chamber under anesthesia (isoflurane) and then hair was removed from the bank flank of the mouse with an electric razor. The mice were then re-anesthetized and placed into a nose cone and injected with buprenorphine in the cephalad region of the surgical site. An incision was created caudally forming a surgical pocket and the titanium disk (coated or uncoated) was inserted. The incision was then closed using surgical staples and mice were monitored until they were conscious.

Histological Analysis

Mice were euthanized using CO₂ inhalation and tissue was harvested (spleen, kidney, liver, intestine, thymus, heart, lung and bone). Tissue samples were fixed, processed, embedded, sectioned (5 microns) and stained with Hematoxylin and Eosin.

Silver Content Analysis

Mice were euthanized using CO₂ inhalation and tissue was harvested (spleen, kidney, liver, intestine, thymus, heart, lung and bone). Each tissue was weighed and 0.5-1 mg of tissue was added to a test tube containing 12 ml of 70% nitric acid. Samples were boiled at 150°C until tissue was liquified and the reaction equalized. Samples were diluted by adding 13 ml of ddH₂O and analyzed by ICP-MS to determine silver content.

AST/ALT Assay

Blood was collected (400 ul) from the facial vein of mice and allowed to coagulated for 2 hours at room temperature. Serum as then used to perform the AST and ALT Elisa assay according to the manufacturer's instructions.

Conclusion

Preliminary data suggests that although there is some evidence of silver accumulation in organs, however, there is no evidence that liver function is compromised. In contrast, histological analysis at later time points shows simple steatosis of the liver at the 1-year time point. Currently, ICP-MS analysis of 1-year tissue is being performed and AST/ALT assays are being repeated for short and long time points. Further histological analysis using autometallography is also being performed to determine the presence of silver in tissue and to complement the ICP-MS analysis.

References

1. Trautner BW, Darouiche RO. Catheter-associated infections: pathogenesis affects prevention. *Arch Intern Med.* 2004;164(8):842-850. doi:10.1001/archinte.164.8.842
2. Biofilm-Related Infections: Bridging the Gap between Clinical Management and Fundamental Aspects of recalcitrance toward Antibiotics. Lebeaux D, Ghigo JM, Beloin C. *Microbiol Mol Biol Rev.* 2014 Sep;78(3):510-43. Review.
3. Silver and Antibiotic, New Facts to an Old Story. Barras F, Aussel L, and Ezraty B. *Antibiotics (Basel).* 2018 Aug 22;7(3). pii: E79. Review.
4. Silver Oxide Coatings with High Silver-Ion Elution Rates and Characterization of Bactericidal Activity. Goderecci S, Kaiser E, Yanakas M, Norris Z, Scaturro J, Oszust R, Medina CD, Waechter F, Heon M, Krchnavek RR, Yu L, Lofland SE, Demarest RM, Caputo GA, and Hettinger JD. *Molecules.* 2017 Sep 7;22(9). pii: E1487.
5. Silver-oxide coating prevents adherence of gram-negative bacteria to titanium biomaterials *in vivo*. Goderecci SS, Cortes E, Manna V, Krchnavek RR, Hettinger JD, Caputo GA, Demarest RM. *In Preparation.*

Acknowledgements

This work was supported by the New Jersey Health Foundation, RowanSOM SMRF and PREP programs and RowanGSBS.